

# *insta***BAG**<sup>®</sup>

Filter blender bag with dehydrated media

## Your diluent ready in minutes!



**NEW**

Available in BPW, Fraser 1/2 and MILK

# interscience

# Ready-to-use media for easy preparation

*instaBAG*<sup>®</sup> is a filter bag which includes a dehydrated pre-dosed media, allowing a 2-in-1 process. Skip the time-consuming media preparation.



1 Place the sample in an *instaBAG*<sup>®</sup>



2 Add sterile and purified water



3 Blend the *instaBAG*<sup>®</sup> with the BagMixer<sup>®</sup>

➔ The sample is ready to be pipetted

## All the advantages of a filter bag

Debris are blocked by the filter during blending

## Long shelf life

Zip pouch with desiccant: optimal conservation after opening



## Just add water

Ready-to-use bag with pre-dosed dehydrated media

Complete dissolution during blending

Non-woven filter identical to the bag filter

PATENTED

SHELF LIFE  
≥ 28 MONTHS

STORAGE  
10°C to 30°C

## Available in BPW, Fraser 1/2 and MILK

*instaBAG*<sup>®</sup> EPT as a diluent for counting quality indicators and *Listeria monocytogenes* and for a non-selective pre-enrichment to detect *Salmonella* and Enterobacteria.


*instaBAG*<sup>®</sup> Fraser 1/2 for primary selective enrichment when testing or counting *Listeria monocytogenes* and *Listeria* spp.

*instaBAG*<sup>®</sup> MILK for a non-selective pre-enrichment for counting quality indicators and detecting *Salmonella* in chocolate products.

## Space saving




100  
*instaBAG*



100  
glass bottles

## No waste



## Less garbage



*instaBAG*  
+ diluent + sample



Media bottle  
+ bag + diluent + sample

## Easy preparation



with *instaBAG*.



with media preparation

## Green

*instaBAG*: no water transported



Diluent pouches



## Long shelf life

*instaBAG*



Media bottle (BPW)



# Technical specifications



	<i>instaBAG</i> <sup>®</sup> BPW 90	<i>instaBAG</i> <sup>®</sup> BPW 225	<i>instaBAG</i> <sup>®</sup> Fraser ½	<i>instaBAG</i> <sup>®</sup> MILK 225	<i>instaBAG</i> <sup>®</sup> MILK 900	<i>instaBAG</i> <sup>®</sup> MILK 2250	<i>instaBAG</i> <sup>®</sup> MILK 3375	
<b>OVERVIEW</b>	<b>Reference</b>	114 090	114 225	114 220	115 022	115 090	115 225	115 337
	<b>Ideal</b>	For pre-enrichment		For selective enrichment	For the analysis of chocolate products			
	<b>Precut for pipetting</b>	✓	✓	✓	✓	✓	-	-
	<b>Marking zone</b>	✓	✓	✓	✓	-	-	-
	<b>Volume index</b>	✓	✓	✓	✓	-	-	-
	<b>Type of bag</b>	Bag with lateral filter						
	<b>Bag composition</b>	Multilayer: multicoated reinforced complex						
	<b>Type of filter</b>	Non-woven						
	<b>Bag volume</b>	400 mL			2000 mL	3500 mL		
	<b>Filter porosity</b>	< 250 µm						
	<b>Type of dehydrated media</b>	Buffered peptone water		Supplemented Fraser ½	Skimmed milk powder			
	<b>Weight of dehydrated media</b>	1.8 g	4.5 g	13 g	21.8 g	87.4 g	218.4 g	327.7 g
	<b>Sample weight</b>	10 g (± 5%)	25 g (± 5%)		100 g (± 5%)	250 g (± 5%)	375 g (± 5%)	
	<b>Volume of water to be added (sterile and purified)</b>	88.2 mL (± 2 %)	220.5 mL (± 2 %)	212 mL (± 2 %)	203.2 mL (± 2 %)	812.6 mL (± 2 %)	2031 mL (± 2 %)	3047 mL (± 2 %)
	<b>Blending time</b>	1 minute		2 minutes	3 minutes			
	<b>Gamma ray treated</b>	10 to 25 kGy, with certificate			25 to 45 kGy, with certificate			
	<b>Shelf-life</b>	45 months		30 months	28 months			
	<b>In compliance with</b>	ISO 11133, ISO 11290-2, ISO 21528, ISO 7218, ISO 6579, ISO 6887 and FDA BAM (Bacteriological Analytical Manual)		ISO 11133, ISO 11290-2, ISO 7218, ISO 6887 and FDA BAM (Bacteriological Analytical Manual)	ISO 7218, ISO 6887 and FDA BAM (Bacteriological Analytical Manual)			
	<b>Approved for food contact</b>	Regulation (EC) No 1935/2004						
	<b>Storage conditions</b>	Temperature from 6 °C to 30 °C (42.8 °F to 86 °F). Close the pouch after opening to avoid humidity		Temperature from 10 °C to 30 °C (50 °F to 86 °F). Close the pouch after opening to avoid humidity	Temperature from 15 °C to 25 °C (59 °F to 77 °F). Close the pouch after opening to avoid humidity			
	<b>Bag dimensions</b>	190 x 300 mm			250 x 380 mm	380 x 600 mm		
	<b>Aluminium pouch of</b>	10 bags			5 bags			
	<b>Box of</b>	100 bags	80 bags	120 bags	40 bags	30 bags	20 bags	
<b>WEIGHT AND DIMENSIONS</b>								

MSDS available on request. Each production batch is checked before shipment (certificate of conformity provided).

**Ask for free samples**

**Certified production**

Products made for INTERSCIENCE by Interlab, an ISO 9001 certified company.



Your local distributor

**interscience**

**PARIS**

Phone: +33 (0)1 34 62 62 61 - Email: info@interscience.com

**FRANKFURT**

Phone: +49 611 7238 7770 - Email: sales.germany@interscience.com

**BOSTON**

Phone: +1 781 937 0007 - Email: sales.usa@intersciencelab.com

**SHANGHAI**

Phone: +86 (0)21-64739390 - Email: sales.china@interscience.cn

**SINGAPORE**

Phone: +65 6977 7232 - Email: sales.asia@interscience.com

**TOKYO**

Phone: +81 3 6712 9715 - Email: sales.japan@interscience.com

[www.interscience.com](http://www.interscience.com)