100 Degrees Celsius Mini Digital Dry Bath ID-SH100-100

The ID-SH100-100 Mini Digital Dry Bath is an extremely compact digital dry bath heat block that provides quick heat up and excellent temperature uniformity and accuracy for sample tubes from 0.2ml up to 50ml. The ID-SH100-100 Mini Digital Dry Bath includes a hinged lid and a block lifter for removing blocks from the unit.

Installation:

Place the ID-SH100-100 Mini Digital Dry Bath on a flat surface away from air vents or exhaust fans from other equipment.

Plug the power pack cord round connector into the back of the Mini Digital Dry Bath. Then connect the power pack plug to a proper grounded electrical outlet.

Open the clear lid and place the aluminum block in the heating well. Load your sample tubes into the block and close the lid. If the sample tubes are too tall, leave the lid open.

Operation:

- Turn the unit on via the ON/OFF switch. The display will illuminate.
- Use the "SET" key to step the cursor through the digits and set each digit to the desired temperature using the up and down arrows. Finish by stepping the cursor back to the left digit.
- The instrument will now quickly heat to the desired temperature. The right side of the display will read "NOT OK" until the desired temperature has been reached. The display will then read "OK".

CAUTION: Hot surfaces, especially on the back, can cause serious injury or burns.

CAUTION: Do not put water or liquids into the well as shock, serious injury and death may occur.

CAUTION: Do not heat flammable or explosive substances as serious injury and death may occur.

Specifications:

Temp. Range:	Ambient +5 to 100°C
Heat Rate:	Under 5 minutes
Temps. Accuracy:	+/- 0.5°C
Temp. Increments:	0.1°C
Temp. Uniformity:	+/- 0.2°C
Block Construction:	High Grade Aluminium
Dimensions:	4.4 x 5.9 x 4 in 11.2 x 15 x 10 cm
Weight:	1.2lbs / 0.54kg
Electrical:	120V to 240V 50-60Hz, 35W
Electrical Optional:	12V adapter for auto outlets
Warranty:	2 years.

Cleaning and Maintenance:

The aluminum blocks can be removed for cleaning by using the included block lifter. Once removed, the blocks can be cleaned with isopropyl alcohol or can be sterilized in an autoclave at 121°C for 20 minutes.

The instrument can be cleaned with a damp cloth.

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