
*F*flame¹¹⁰ Safety Bunsen

NEW
Safety Bunsen

Flaming from any angle

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WLD-TEC GmbH

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Bracket sidewise (r/l)
stainless steel

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


Bracket horizontal
stainless steel

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* optional

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Flame¹¹⁰ – the extraordinary precision burner

Are you searching for a solution to heat or flame any item under any angle - even upside down?

Here's what you need:

Our new precision burner Flame¹¹⁰!

Multi-talent

The new Flame¹¹⁰ is ideally suited for any kind of flaming process in laboratory or industry.


Its long precision flame can heat your dental instruments, can sterilize your microbiological tools or flame deburr resp. flame polish your moulded plastic or rubber parts in any type of industrial production. Always safe and perfect!

Simple and safe handling

Simply push the button: Flame will start or stop instantly! You may also control it by foot pedal* or external IR motion sensor*. Safety is built-in and always active (overheat protection and automatic re-ignition).

Compact and efficient

Stainless steel, efficient cleaning, small footprint, exclusive variability and convincing economics are vital for any type of laboratory and industrial workstation. Flame¹¹⁰ is offering a perfect combination of features:

-  variable flame from any angle – even upside down
-  safety in any operation mode
-  very easy cleaning of removable burner head
-  operation with liquid gas from cylinders or cartridges (with optional adapters)
-  very economic consumption of gas

- (1) **Fußpedal / Foot pedal**
Edelstahl / stainless steel
Art.-Nr. / Art.-No. 6.000.402
- (2) **Funk-Fußpedalset mit RF-Stick**
(nur in EU Staaten lieferbar)
Radio controlled foot pedal with RF-Stick
(available in EU countries only)
Art.-Nr. / Art.-No. 6.000.404-RF
- (3) **Fußpedal Mini / Foot pedal Mini**
Kunststoff / plastic
Art.-Nr. / Art.-No. 6.000.403

Externer DoubleClick IR-Bewegungssensor
zum berührungslosen
Flammenstart, Edelstahl,
Kabellänge 0,4 m
External DoubleClick IR-motion sensor
touch-free flame activation,
stainless steel, cable length 0.4 m
Art.-Nr. / Art.-No. 6.000.406

Flame¹¹⁰
Safety Bunsen

Programs

Button function	✓
Pedal Standard (food pedal optional)	✓
Gas inlet filter	✓
Removeable burner head	✓
Housing, gas connection & controls made of stainless steel	✓
DVGW approved	✓

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Safety Bunsen

Art.-No. 2.110.000-L

Nozzle for propane / butane gas
Removable special burner head
Bracket with 2 wing nuts for
upside down mounting
Wrench 17 mm
Switching power supply (worldwide,
level 6)
Instruction manual

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Technical Data:

Programs

Button: Start-Stop with timer, 60 min
flame during pressed foot pedal

Foot pedal (optional):

Safety features

with gas safety cut off:

re-ignition 1 x per second
overheat protection
pressure monitoring at approx. 3.5 to 4 bar

Gas supply and consumption

Gas connection: 1/4" left hand thread with gas filter
Gas types: I3P 1.5 bar liquid gas
Connected load: 38 g/h
Cartridge operation (type - time): C 206 - 4,5 h / CV 300 - 6 h / CV 470 - 11 h

Flame characteristics

Max. flame temperature: approx. 1200 °C
Flame measurements (h x Ø): max. 120 x 10 mm
Temperature threshold level: 500 W

Electrical

Power consumption: 2 VA (Stand-by max 0,1 VA)
Switching power supply: 100-240V / 50/60Hz / max. 0,3A; 9V DC / 1,3A
Level 6 (Stand-by Verbrauch max. 0,1 W)

Mechanical

Burner head: removable

Dimensions (h x Ø) / with burner head: 75 x 88 mm / 135 x 88 mm
Weight: 580 g

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Safety adapter for gas cartridges
with pressure regulator and DVGW-tubing (0.5 m)

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
Safety medium pressure regulator - 1.5 bar
for gas cylinders of 5, 11 or 33 kg

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DVGW Safety tubing
for threaded connection 3/8" L«-» 1/4" L (1.5 m)

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