U.S. English

Product Number:
902086G

Reveal[®] 3-D
for Almond

Reveal[®] 3-D for Almond

Product Number: 902086G

Introduction and Intended Use

Reveal[®] 3-D for Almond is uniquely designed with three lines of detection and can be used virtually anywhere to screen environmental swabs, clean in place (CIP) rinses, food products and ingredients. For the presence of significant levels of almonds. The test's 3-D technology ensures greater reliability with screening than ever before.

Reveal 3-D for Almond is intended for use only in an industrial food manufacturing/preparation context or for food labelling enforcement testing. Because of the problems of adequately sampling and extracting almond, it is not suitable for the testing of foods to be consumed in the home or in a restaurant by allergic individuals. The test detects significant (low parts per million, or ppm) levels of almond content in CIP rinses, environmental swabs, and foods.

Detection Limit

The test uses highly specific antibodies to detect almond protein. Almond can be detected from various surfaces using the provided environmental swabs, rinses and liquid products.

Utilizing the environmental swabs supplied, levels of 1 µg/ 100 cm² of almond on surfaces can be detected.

When analyzing CIP rinses and foods, almond residues are detectable as low as 1.1 ppm almond protein (5 ppm total almond). The presence of cleaners and sanitizers in rinses, as well the type of food matrix being tested, can affect limit of detection.

Cross-reactivity

The test was used to analyze a panel of potentially cross reacting commodities including nuts, seeds, legumes, and grains. This test cross-reacts with common *Prunus* genus (i.e. apricots, cherries, nectarines, peaches and plums). Cross-reactivity has been identified with the seed (stone); however, the flesh (seedless fruit) does not cross react. Further information is detailed in the Reveal 3-D Almond customer validation report, which is available on request.

Test Performance

Reveal 3-D for Almond has undergone rigorous validation to evaluate the specificity, sensitivity, robustness and intra- and inter-batch variability of the test method on CIP rinses, environmental swabs and foods. A customer validation report is available on request.

Sample Compatibility

Reveal 3-D for Almond is designed to detect almond on environmental surfaces, CIP rinses, and foods. Although every effort has been made to validate as many variables as possible, there may be some sample types that are not suitable for testing.

Users should perform in-house matrix-specific spike recovery validation work in conjunction with a validated laboratory assay (e.g., Veratox[®]) to help confirm Reveal 3-D for Almond results. This procedure will highlight any problematic matrices encountered.

For certain food ingredients or more complex food matrices, it is recommended that a more broadly inclusive test, such as Veratox for Almond, be used. Existing commodity validations should be revalidated when suppliers or the manufacturing process has been changed. Please contact a NEOGEN representative or distributor for additional details.

Reagents and Materials Provided

The Reveal 3-D Almond test pack contains the following:

1. An instruction leaflet
2. 1 foil pouch, containing 10 yellow Reveal 3-D for Almond devices
3. 10 sachets containing Type 12 Extraction Buffer
4. 10 sample tubes and caps
5. 10 individually packaged sterile swabs with break off tips
6. 1 bottle of swab wetting solution

Materials Required for Food Extraction but not Provided

1. 3-D Food Buffer (NEOGEN item 8504)
2. 50 mL graduated cylinder tubes for extraction (NEOGEN item 9381)

Sampling Technique

1. CIP rinse water and liquid samples
As only a small amount of material (0.25 mL) is required for the Reveal 3-D for Almond test, it is important to test a representative portion of liquid.
2. Environmental swabbing
The swabs supplied are intended to be used for the collection of environmental samples from which the presence of almond can be tested. This method can be used to validate the adequacy of cleaning and/or to identify problem areas (e.g., unwanted build-up of almond in processing equipment).
3. Food samples
As only a small sample of material is required for Reveal 3-D for Almond, it is important to test a representative portion of the sample.

Precautions

1. For testing in an industrial food manufacturing/ preparation, or labeling enforcement context only.
2. Do not use any part of the test beyond the expiry date.
3. Do not open the foil pouch until just before use.
4. Always store the kit between 2–8°C (35–46°F). Avoid freezing.
5. Ensure the foil bag is tightly sealed after removal of a device.

Test Storage / Sample Stability


Store the Reveal 3-D for Almond kit between 2–8°C (35–46°F). The kit should be used within the expiry date stated on the outer label. Extracted samples should be used within 3 hours of extraction.

Limitations


A negative test from a surface or rinse cannot exclude the possibility that the food contains almond since it may be distributed unevenly on the surface or sample and may be below the detection limit of the test with that specific sample.

Reveal 3-D for Almond is qualitative and should only be used as a preliminary screen for the presence of almond. The validity of results obtained with the test should preferably be viewed in conjunction with data from a quantitative assay such as Veratox for Almond allergen.

Validate detection by testing a positive control from the source allergen or ingredient that is available in your manufacturing environment and that may pose an allergen contamination risk to ensure the test can detect the allergen of concern. Contact NEOGEN Technical Services for additional information.

1


2


Rinse Sampling

Remove the following and allow to equilibrate at room temperature before use (20–30 minutes out of refrigerator):

1. 1 sample tube
 2. 1 Type 12 Extraction Buffer sachet
 3. 1 Reveal 3-D for Almond device (in foil pouch)
1. Carefully tear/cut off the top of the clear Type 12 Extraction Buffer sachet and add the entire contents to the sample tube.
 2. Add 0.25 mL (250 ul) of sample to the sample tube.
 3. Secure the white cap and shake for 1 minute.

1

2


Swab Sampling

Remove the following and allow to equilibrate at room temperature before use (20–30 minutes out of refrigerator):


1. 1 bottle swab wetting solution
 2. 1 sample tube
 3. 1 Type 12 Extraction Buffer sachet
 4. 1 Reveal 3-D for Almond device (in foil pouch)
 5. 1 sterile swab
1. Carefully tear/cut off the top of the Type 12 Extraction Buffer sachet and add the entire contents to the sample tube.
 2. Estimate a swabbing area of approximately 10 x 10 cm. Alternatively, use the swab to collect samples of contamination from problem areas (e.g., of processing equipment, filler heads).
 3. Gather the sample with the swab, using one of the following methods:

For dry surfaces: Remove a sterile swab from the packaging and wet with two drops of swab wetting solution. Swab a 10 x 10 cm area by using a crosshatch technique revolving the swab on the surface. Repeat this swabbing procedure using movements at right angles to those used in the first swabbing.

For wet surfaces: Remove a sterile swab from the packaging and swab a 10 x 10 cm area by using a crosshatch technique, revolving the swab on the surface. Repeat this swabbing procedure using movements at right angles to those used in the first swabbing. Do not moisten swab prior to use.


3



4

4. Return the swab to the Extraction Buffer in the sample tube and carefully break off the moistened end at the pre-scored mark so it remains in the tube.
5. Secure the cap of the sample tube, taking care to ensure that the stem does not prevent the tube from being properly sealed. Shake for 1 minute.

5




NEOGEN.com

Swab and Rinse Testing

1. Remove the lid and fill it with the liquid from the tube. Any froth should remain in the tube.
2. Dip the head of the Reveal 3-D device into the liquid in the lid. Ensure that the cavity is saturated with the liquid.
3. Leave the cavity saturated until you see the liquid running in the test window.
4. Place device on a flat surface and allow test to develop for 5 minutes.

Note: It is essential to place the device flat on a level surface as soon as the liquid has entered the test window to stimulate flow through the device. Additionally, the devices are pre-stripped with pale green loading dye in positions T (test), O (overload), and C (control). The loading dye is removed from the test window as the sample flows through the device.





Food Extraction and Testing

The extraction of food samples requires the use of NEOGEN's 3-D Food Buffer (NEOGEN item 8504).

From the test kit, remove the appropriate number of the following and allow to equilibrate at room temperature before use (20–30 minutes out of refrigerator):

1. 3-D Food Buffer
 2. Reveal 3-D for Almond device (in foil pouch)
1. Add 20 mL of 3-D Food Buffer into a new 50 mL tube (18 mL for liquid samples) or a similar container.
 2. Weigh out 2 g (2 mL for liquid products) of sample and add to tube with buffer.
 3. Secure the caps and invert to suspend sample. Shake vigorously or vortex for 1 minute.
 4. Tilt the tube at approximately a 45 degree angle. Dip the head of the Reveal 3-D device into the liquid. Ensure the cavity is saturated with the liquid, and that no liquid gets on the test window.
 5. Leave the cavity saturated until liquid is seen running in the test window
 6. Place device on a flat surface and allow test to develop for 5 minutes.

Note: It is essential to place the device flat on a level surface as soon as the liquid has entered the test window to stimulate flow through the device. Additionally, the devices are pre-stripped with pale green loading dye in positions T (test), O (overload), and C (control). The loading dye is removed from the test window as the sample flows through the device.

Reading Results

Read the result 5 minutes after dipping. Observations after 6 minutes may be inaccurate due to overdevelopment of the device.

1. Negative result
No line at position T (test): Level of almond is below the limit of detection. (See Limitations section.)
2. Positive result
Any intensity of the red line at position T (test): Level of almond above detection limit.
3. High positive results
No line is visible at position O (overload) and a red line may be faintly visible or absent at position T: Sample is overloaded with almond.
4. Invalid results
If no line appears at position C (control), then the test may be invalid.

Customer Service

NEOGEN Customer and Technical Services can be contacted through [NEOGEN.com](https://www.neogen.com) and product training is available by request.

SDS Information Available

Safety data sheets are available for all test kits at [NEOGEN.com](https://www.neogen.com) or by calling 800.234.5333 or 517.372.9200.

Terms and Conditions

NEOGEN's full terms and conditions are available [online](#).

Warranty

NEOGEN makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, NEOGEN will provide a replacement of the product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. NEOGEN shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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