U.S. English

Product Number:
901041L

Reveal 3-D
for Peanut

Reveal[®] 3-D for Peanut

Product Number: 901041L

Introduction and Intended Use

Reveal[®] 3-D for Peanut is uniquely designed with three lines of detection and can be used virtually anywhere to screen environmental swabs, clean-in-place (CIP) rinses, food products and ingredients for the presence of significant levels of peanut. The test's 3-D technology ensures greater reliability with screening than ever before.

Reveal 3-D for Peanut is intended to be used only in an industrial food manufacturing or preparation context. Because of the problems of adequately sampling and extracting peanut, it is not suitable for the testing of foods to be consumed in the home or in a restaurant by allergic individuals. The test detects significant (low parts per million, or ppm) levels of peanut content in CIP rinses, environmental swabs, and foods.

Limit of Detection

The test uses highly specific antibodies to detect allergenic proteins from peanuts. Residues of peanuts can be detected from various surfaces using the provided environmental swabs or from rinses and foods using the below testing protocols.

Using the environmental swabs supplied, levels below 3–4 µg/100 cm² of peanut on surfaces can be detected.

When analyzing CIP rinses and foods, peanut residues are detectable as low as 1.3 ppm peanut protein (5 ppm total peanut). The presence of cleaners and sanitizers in rinses, as well the type of food matrix being tested, can affect limit of detection.

Test Performance

Reveal 3-D for Peanut has undergone rigorous validation to evaluate the specificity, sensitivity, robustness and intra- and inter-batch variability of the test method on CIP rinses, environmental swabs, and foods. A customer validation report is available on request. In an AOAC Research Institute Performance Tested Method License # 111901 study, the Reveal 3-D Peanut test was found to be an effective procedure for qualitative determination of the presence of peanut residues in environmental swabs and CIP rinses.

Sample Compatibility

Reveal 3-D for Peanut is designed to detect peanut on environmental surfaces and in CIP rinses, and foods. Although every effort has been made to validate as many variables as possible, there may be some sample types that are not suitable for testing. The test is best used in conjunction with a validated laboratory assay (e.g., Veratox[®] for Peanut allergen) to help confirm Reveal 3-D for Peanut results.

Users should perform in-house matrix-specific spike and recovery validation work in conjunction with a validated laboratory assay (e.g. Veratox) to help confirm Reveal 3-D for Peanut results. This process will highlight any problematic matrices encountered.

Existing commodity validations should be revalidated when suppliers or the manufacturing process has been changed. Please contact a NEOGEN representative or distributor for additional details.

Materials Provided

The Reveal 3-D for Peanut kit contains the following:

1. An instruction leaflet
2. 1 foil pouch, containing 10 orange Reveal 3-D for Peanut test devices
3. 10 bottles or sachets containing Type 9 Extraction Buffer
4. 10 sample tubes and caps
5. 10 individually packaged, sterile swabs with break-off tips
6. 1 bottle of swab wetting solution

Materials Required For Food Extraction But Not Provided

1. 3-D Food Buffer (NEOGEN item 8504)
2. 50 mL graduated cylinder tubes for extraction (NEOGEN item 9381)

Sampling Technique

1. CIP rinse water
As only a small amount of material is required for Reveal 3-D for Peanut, it is important to test a representative portion of liquid.
2. Environmental swabbing
The swabs supplied are intended to be used for the collection of environmental samples from which the presence of peanut can be tested. This method can be used to validate the adequacy of cleaning and/or to identify problem areas (e.g., unwanted buildup of peanut in processing equipment, filler heads).
3. Food samples
As only a small sample of material is required for Reveal 3-D for Peanut, it is important to test a representative portion of the sample.

Precautions

1. For environmental testing in an industrial food manufacturing/preparation, or labeling enforcement context only.
2. Do not use any part of the test beyond the expiry date.
3. Do not open the foil pouch until just before use.
4. Ensure the foil bag is tightly sealed after removal of a device.
5. Always store the kit between 2–8°C (35–46°F). Avoid freezing.
6. Bring kit to room temperature 18–30°C (64–86°F) prior to use.

Test Storage/Sample Stability

Store the Reveal 3-D for Peanut kit between 2–8°C (35–46°F), and use within the expiry date stated on the outer label. Extracted samples should be used within 3 hours of extraction.

Limitations


A negative test from a surface, rinse or food cannot exclude the possibility that the sample contains peanut since it may be distributed unevenly on the surface or rinse and may be below the detection limit of the test.

Reveal 3-D for Peanut is qualitative and should only be used as a preliminary screen for the presence of peanut. The validity of results obtained with the test should preferably be viewed in conjunction with data from a quantitative assay, such as Veratox for Peanut allergen.


Samples tested containing >10% walnut may interfere with the assay and the method may not be able to detect peanut residues. Due to the nature of the samples, environmental swabs and rinses typically will have concentrations of walnut well below this threshold. Please contact a NEOGEN representative for further details.

Validate detection by testing a positive control from the source of the allergen or ingredient that is available in your manufacturing environment and that may pose an allergenic contamination risk to ensure the test can detect the allergen of concern. Contact NEOGEN technical services for additional information.

NEOGEN.com

1


3


Rinse Sampling

For each sample to be tested, remove the following and allow to equilibrate at room temperature before use (20–30 minutes out of refrigerator):


1. 1 sample tube and cap
 2. 1 Type 9 Extraction Buffer
 3. 1 Reveal 3-D for Peanut device (in foil pouch)
1. Carefully tear/cut or uncap the Type 9 Extraction Buffer and add the entire contents to the sample tube.
 2. Add 0.25 mL (250 ul) of sample to the sample tube.
 3. Secure the white cap and shake for 1 minute.


1


2

3

4

5

Swab Sampling


For each sample to be tested, remove the following and allow to equilibrate at room temperature before use (20–30 minutes out of refrigerator):

1. 1 bottle swab wetting solution
 2. 1 sample tube
 3. 1 Type 9 Extraction Buffer
 4. 1 Reveal 3-D Peanut Test device (in foil pouch)
 5. 1 sterile swab
1. Carefully tear/cut the Type 9 Extraction Buffer and add the entire contents to the sample tube.
 2. Estimate a swabbing area of approximately 10 x 10 cm.
Alternatively, use the swab to collect samples of contamination from problem areas (e.g., of processing equipment, filler heads).
 3. Gather the sample with the swab, using one of the following methods:

For dry surfaces: Remove a sterile swab from the packaging and wet with 2 drops of swab wetting material. Swab a 10 x 10 cm area by using a crosshatch technique, revolving the swab on the surface. Repeat this swabbing procedure using movements at right angles to those used in the first swabbing.

For wet surfaces: Remove a sterile swab from the packaging and swab a 10 x 10 cm area by using a crosshatch technique, revolving the swab on the surface. Repeat this swabbing procedure using movements at right angles to those used in the first swabbing. Do not moisten swab prior to use.



4. Return the swab to the Extraction Buffer in the sample tube and carefully break off the moistened end at the pre-scored mark so that it remains in the tube.
5. Secure the cap of the sample tube, taking care to ensure that the stem does not prevent the tube from being properly sealed. Shake for 1 minute.

Swab and Rinse Testing

1. Remove the lid and fill it with the liquid from the tube. Any froth should remain in the tube.
2. Dip the head of the Reveal 3-D device into the liquid in the lid. Ensure that the cavity is saturated with the liquid.
3. Leave the cavity saturated until liquid is seen running in the test window.
4. Place device on a flat surface and allow test to develop for 5 minutes.

Note: It is essential to place the device flat on a level surface as soon as the liquid has entered the test window to stimulate flow through the device. Additionally, the devices are pre-stripped with pale green loading dye in positions T (test), O (overload), and C (control). The loading dye assists with quality and manufacturing checks and does not impact test performance. The loading dye is removed from the test window as the sample flows through the device.

Food Extraction and Testing

The extraction of food samples requires the use of NEOGEN's 3-D Food Buffer (NEOGEN item 8504).



From the test kit, remove the appropriate number of the following and allow to equilibrate at room temperature before use (20–30 minutes out of refrigerator):





3-D Food Buffer

Reveal 3-D for Peanut device (in foil pouch)

1. Add 20 mL of 3-D Food Buffer into a new 50 mL tube (18 mL for liquid samples) or a similar container.
2. Weigh out 2 g (mL) of sample and add to tube with buffer.
3. Secure the caps and invert to suspend sample. Shake vigorously or vortex for 1 minute.
4. Tilt the tube at approximately a 45 degree angle. Dip the head of the Reveal 3-D device into the liquid. Ensure the cavity is saturated with the liquid, and that no liquid gets on the test window.
5. Leave the cavity saturated until liquid is seen running in the test window
6. Place device on a flat surface and allow test to develop for 5 minutes.

Note: It is essential to place the device flat on a level surface as soon as the liquid has entered the test window to stimulate flow through the device. Additionally, the devices are pre-stripped with pale green loading dye in positions T (test), O (overload), and C (control). The loading dye assists with quality and manufacturing checks and does not impact test performance. The loading dye is removed from the test window as the sample flows through the device.

Reading Results

Read the result 5 minutes after dipping. Observations after 6 minutes may be inaccurate due to overdevelopment of the device.

1. Negative result
No line at position T (test): Level of peanut is below the detection limit. (See Limitations section).
2. Positive result
Any intensity of a red line at position T (test): Level of peanut above detection limit.
3. High positive results
No line is visible at position O (overload) and a red line is faintly visible or absent at position T: Sample is overloaded with peanut.
4. Invalid results
If no red line appears at position C (control), then the test may be invalid.

Customer Service

NEOGEN Customer and Technical Services can be contacted through NEOGEN.com and product training is available by request.

SDS Information Available

Safety data sheets are available for all test kits at NEOGEN.com or by calling 800.234.5333 or 517.372.9200.

Terms and Conditions

NEOGEN's full terms and conditions are available [online](#).

Warranty

NEOGEN makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, NEOGEN will provide a replacement of the product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. NEOGEN shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product. fitness of the product for any purpose. NEOGEN shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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