*Read instructions carefully before starting test*

# **Reveal<sup>®</sup>**

## **for CP4 (Roundup Ready<sup>®</sup>)**

### **CP4**

Monsanto's Roundup Ready<sup>®</sup> crops are plants that have been genetically improved by the insertion of the CP4 EPSPS coding sequence. The expression of the CP4 protein makes the crops tolerant to Roundup<sup>®</sup> herbicides, allowing an over-the-top application of the herbicide to control weeds without endangering the safety of the crop.

### **INTENDED USE**

Reveal for CP4 is an immunochromatographic lateral flow assay used for the qualitative analysis of CP4 in soybeans and corn. The test can detect as little as 1 CP4 enhanced soybean among 1,000 beans (0.1%), and 1 NK603 Roundup Ready kernel of corn among 800 kernels (0.125%). Earlier Roundup Ready corn varieties, such as GA21, will not be detected. This test is not valid for testing wheat or wheat products.

### **INTENDED USER**

The test kit is designed for use by quality control personnel and others familiar with crops possibly containing CP4. Since technique is very important, operators should be trained by a Neogen representative or someone who has successfully completed Neogen training.

### **ASSAY PRINCIPLES**

Reveal for CP4 is a single-step lateral flow immunochromatographic assay. The water extract is wicked through a reagent zone, which contains antibodies specific for CP4 conjugated to colored particles. If CP4 is present, it will be captured by the conjugated antibodies. The CP4-antibody-particle complex is then wicked onto a membrane which contains a zone of antibody specific for CP4. This zone captures the complex allowing the particles to concentrate and form a visible line. If no CP4 is present, no line will form. The membrane also contains a control zone where an immune complex present in the reagent zone is captured by an antibody, forming a visible line. The control line will always form regardless of the presence of CP4, ensuring the strip is working properly.

## STORAGE REQUIREMENTS

Store kit components at room temperature (18–30°C/64–86°F) to assure full shelf life.

## MATERIALS PROVIDED

1. 50 CP4 test strips
2. 50 transfer pipettes
3. 50 sample tubes

## MATERIALS REQUIRED BUT NOT PROVIDED

1. 32 oz Mason-type (glass or plastic) jars
2. Blender blade assembly for use with 32 oz jars
3. Scale capable of weighing 300 grams (Neogen item 9427)
4. Distilled, deionized or tap water
5. Timer (Neogen item 9426)
6. Sample tube rack (Neogen item 9466)
7. Waterproof marker
8. Graduated cylinder
9. Centrifuge, mini (Neogen item 9330)
10. Tubes, 1.5 mL micro-centrifuge tubes – box of 500 (Neogen item 9372)

## PRECAUTIONS

1. Store test kit at room temperature (18–30°C/64–86°F) when not in use.
2. Do not use kit components beyond expiration date.
3. Do not freeze test kits.
4. Treat all used liquids, including sample extract, and labware as if containing CP4.
5. To avoid cross-contamination, use clean labware for each sample, and thoroughly wash all labware between samples.

## SOYBEAN SAMPLE PREPARATION AND EXTRACTION

The sample to be tested should be collected according to accepted sampling techniques. Extracted samples should be tested immediately.

1. Determine the average bean weight of the soybeans to be tested by counting out 100 beans, weighing all 100 together in one container to the nearest 0.01 gram, and dividing the total by 100.

**NOTE:** The result in step 1 should approximately correlate with the size classifications for soybeans established by the Association of Official Seed Analysts. The classifications are: Very Petite, 0.076 g/seed; Small, 0.095 g/seed; Medium, 0.223 g/seed; and Large, 0.293 g/seed.

2. To achieve 0.1% sensitivity with this test, multiply the average bean weight determined in step 1 by 1,000 to determine the sample size. (Example: An average bean weight of 0.095 X 1,000 = a sample size of 95 grams.)
3. Weigh out the soybean sample size determined in step 2 and deposit into a 32 oz jar. Attach blender blade assembly and blend on high speed for 30-45 seconds.
4. Remove the blender blade and add 4 times the sample size of water to the jar. (Example: 4 X 95 grams of sample = 380 mL of water.) Cap and shake for 30-45 seconds to ensure all the blended soybean has been wetted.
5. Let soybean and water mixture settle, and use liquid from the top of the jar as the test sample.

**NOTE:** It is highly recommended to centrifuge soy samples by pipetting 1 mL of sample into a 1.5 mL micro-centrifuge tube, and centrifuge for 1 minute.

## CORN SAMPLE PREPARATION AND EXTRACTION

The sample to be tested should be collected according to accepted sampling techniques. Extracted samples should be tested immediately.

1. Weigh 250 g  $\pm$  10 g of whole corn sample to be tested and deposit into a 32 oz jar. Attach blender blade assembly and blend on high speed for 30–45 seconds.
2. Remove blender blade and add 400  $\pm$  20 mL of water to the jar. Cap and shake for 30–45 seconds to ensure all corn has been wetted.
3. Let corn and water mixture settle, and use liquid from the top of the jar as the test sample.

## TEST PROCEDURE

1. Remove the appropriate number of sample tubes and place in tube rack.  
**IMPORTANT:** You must use the supplied sample tubes for the test procedure for proper sample flow.
2. Using a new transfer pipette for each sample, transfer approximately 0.5 mL of extract to a sample tube.
3. Place a new CP4 strip with the sample end down into a sample tube.
4. Allow the strip to develop and read results.

## INTERPRETATION OF RESULTS

**Positive results:** If a line forms in the test zone and another line forms in the control zone within 10 minutes (2 lines total), the sample is positive for CP4 at  $\geq$  0.1% for soybeans, or  $\geq$  0.125% for corn. A positive sample may be determined as soon as 2 lines are visible on the strip. If there is no line in the control zone, the test is invalid and the sample should be re-tested with another strip.

**Negative results:** If after the full 10 minutes there is no visible line in the test zone, but a visible line in the control zone, the sample is negative for CP4. If there is no line in the control zone, the test is invalid and the sample should be retested with another strip.

## NOTES:

- Positive results can be recorded as soon as 2 lines are visible, within 1–10 minutes of development time.
- Negative results must be observed after 10 minutes. Observations after 15 minutes may be inaccurate due to over-development of the device.

*\*Contact Neogen or visit [www.neogen.com](http://www.neogen.com) for more information on these products.*

**CUSTOMER SERVICE**

Neogen Customer Assistance and Technical Services can be reached by using the contact information on the back of this booklet. Training on this product, and all Neogen test kits, is available.

**SDS INFORMATION AVAILABLE**

Safety data sheets (SDS) are available for this test kit, and all of Neogen's Food Safety test kits, on Neogen's Web site at [www.neogen.com](http://www.neogen.com), or by calling Neogen at 800/234-5333 or 517/372-9200.

**TERMS AND CONDITIONS**

For Neogen's full terms and conditions, please visit [www.neogen.com/Corporate/termsconditions.html](http://www.neogen.com/Corporate/termsconditions.html).

**WARRANTY**

Neogen Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, Neogen will provide a replacement product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

## TESTING KITS AVAILABLE FROM NEOGEN

### Natural toxins

- Aflatoxin, DON, ochratoxin, zearalenone, T-2/HT-2 toxins, fumonisin, histamine

### Foodborne bacteria

- *E. coli* O157:H7, *Salmonella*, *Listeria*, *Listeria monocytogenes*, *Campylobacter*, *Staphylococcus aureus*, *Salmonella enteritidis*

### Sanitation

- ATP, yeast and mold, total plate count, generic *E. coli* and total coliforms, protein residues

### Food allergens

- Almonds, crustaceans, eggs, gliadin, hazelnut, milk, mustard, peanuts, sesame, soy, walnuts, multi-treenut

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#### India

**Neogen Food and Animal Security**  
+91 484 2306598, 2301582  
info@neogenindia.com  
www.neogenindia.com

©Neogen Corporation, 2017. Neogen, AccuScan and Reveal are registered trademarks of Neogen Corporation. All other trademarks are the property of their respective companies. Patent: <http://www.neogen.com/Corporate/patents.html>.

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