
CONFIRMATION AGARS FOR *LEGIONELLA* (BCYE α AND BCYE α WITHOUT CYSTEINE)

CONFIRMATION FOR ENUMERATION AND ISOLATION OF *LEGIONELLA*

BM07208
BM07308

1 INTENDED USE

The two confirmation agars for *Legionella*, BCYE α and BCYE α without cysteine, are destined for the confirmation of colonies previously isolated on GVPC agar. Collectively, these media are used for the enumeration and isolation of *Legionella* species in water and other samples susceptible of harboring the bacteria.

2 HISTORY

- In 1977, MacDade *et al.* were the first to isolate the agent responsible for Legionnaire's Disease, a bacteria now known as *Legionella pneumophila*. After this discovery, numerous occurrences of *Legionella* isolation of *Legionella* were reported in fresh water environments such as water distribution systems, air conditioning, cooling towers, and spas. 48 species of *Legionella* are currently known.
- In 1978, Weaver succeeded in cultivating *Legionella* on Mueller-Hinton chocolate agar. Feeley *et al.*, deduced that cysteine and ferric pyrophosphate could replace the vitamin and hemoglobin supplements found in the Mueller Hinton chocolate agar. Their work led to the formulation of a medium dubbed F-G agar. They determined as well that an atmosphere enriched at 2.5% CO₂ was necessary for *Legionella* culture.
- In 1979, Feeley *et al.* modified the F-G medium by replacing acid hydrolysate of casein by yeast extract, and adding activated charcoal while eliminating starch. The resulting CYE media allowed better growth of *Legionella*.
- In 1980, Pasculle *et al.* supplemented the CYE medium with ACES buffer. They demonstrated that this new medium, designated BCYE, offered a better recovery of *Legionella* and could be incubated aerobically.
- End in 1981, Edelstein increased the sensitivity of the medium by adding α -cetoglutarate (BCYE α medium). Colonies are considered to be *Legionella* if they develop on BCYE α medium, but demonstrate no growth on BCYE α without cysteine.

3 PRINCIPLES

- Yeast extract constitutes a primary nutrient leading to *Legionella* growth.
- Activated charcoal decomposes hydrogen peroxide (toxic metabolic by-product), captures the carbon dioxide and modifies the surface tension.
- The ACES/KOH buffer maintains the pH and permits aerobic incubation.
- Cysteine and ferric pyrophosphate represents indispensable nutritive elements for the growth of *Legionella*.
- α -cetoglutarate is a growth activator for *Legionella*.

4 INSTRUCTIONS FOR USE

- Select at least three characteristic *Legionella* colonies from each plate of GVPC agar (**BM071**).
- Re-streak each colony onto a plate of BCYE α agar without cysteine (**BM073**) and a plate of BCYE α agar (**BM072**).
- Incubate at 36 \pm 1°C for at least 2 days.
- Consider as positive for *Legionella* all colonies that develop on BCYE α agar but present no growth on BCYE α agar without cysteine.
- Proceed with serological identification of the *Legionella* species.

5 RESULTS

The present of colonies on BCYE α Agar with cysteine (**BM072**) and the absence of colonies on BCYE α Agar without cysteine (**BM073**) proves the present of *Legionella*.

Colonies of *Legionella* spp. present a white to gray coloration on BCYE α Agar with cysteine (**BM072**). They can also have blue, pink, purple, maroon, greenish-yellow or dark red pigmentation that fades, becoming whiter and filamentous with age. Their surface is smooth with precise edges. Some strains may give a ground glass or "fried egg" aspect when observed through a binocular scope, while others may present a brilliant white fluorescence under a UV light. (see APPENDIX 1 : PHOTO SUPPORT)

6 TYPICAL COMPOSITIONS

Typical compositions can be adjusted to obtain optimal performance.

• BCYE α Agar with cysteine

For 1 liter of media :

- Yeast extract	10,0 g
- Activated charcoal	2,0 g
- α -cetoglutarate, monopotassium salt	1,0 g
- ACES (2-[2-amino-2-oxoethyl)-amino] ethanesulfonic acid)	10,0 g
- Potassium hydroxide	2,8 g
- Ferric pyrophosphate	0,25 g
- Bacteriological agar	12,0 g
- L-cysteine, hydrochloride	0,4 g

pH of the ready-to-use media at 25°C : 6.9 \pm 0.1

• BCYE α agar without cysteine

For 1 liter of media :

- Yeast extract	10,0 g
- Activated charcoal	2,0 g
- α -cetoglutarate, monopotassium salt	1,0 g
- ACES (2-[2-amino-2-oxoethyl)-amino] ethanesulfonic acid)	10,0 g
- Potassium hydroxide	2,8 g
- Ferric pyrophosphate	0,25 g
- Bacteriological agar	12,0 g

pH of the ready-to-use media at 25°C : 6.9 \pm 0.2

7 QUALITY CONTROL

BCYE α agar (BM072) / BCYE α agar without cysteine (BM073)

- Prepared media in plates : black agar, with visible particles of activated charcoal.
- Typical cultural response after 72 hours incubation at 36°C :

Microorganisms		Growth	
		BCYE α agar	BCYE α agar without cysteine
<i>Legionella pneumophila</i>	WDCM 00107	positive	negative
<i>Escherichia coli</i>	WDCM 00179	positive	positive

8 STORAGE

Pre-poured media (complete) in Petri dishes BM072 and BM073:

Store between 2 - 8°C, shielded from light.
The expiration date is indicated on the label

9 PRESENTATION

Pre-poured media in Petri plates (Ø 90 mm):

BCYE α agar – 20 plates	BM07208
BCYE α agar without cysteine – 20 plates	BM07308

10 BIBLIOGRAPHY

McDade, J.E., Shepard, C.C., Fraser, D.W. et al. 1977. Legionnaires' disease. Isolation of a bacterium and demonstration of its role in other respiratory disease. N. Engl. J. Med., 297 : 1197-1203.

Feeley, J.C., Gorman, G.W., Weaver, R.E., Mackel, D.C., and Smith, H.W. 1978. Primary isolation media for the Legionnaires' disease bacterium. J. Clin. Microbiol., 8 : 320-325.

Weaver, R.E., and Feeley, J.C. 1979. Cultural and biochemical characterisation of legionnaires disease bacterium. 'Legionnaires', the bacterium and methodology. Center for Disease Control, Atlanta GA. p 19-25.

Feeley, J.C., Gibson, R.J., Gorman, G.W., Langford, N.C., Rasheed, J.K., Mackel, D.C. and Baine, W.B. 1979. Charcoal-yeast extract agar : primary isolation medium for *Legionella pneumophila*. J. Clin. Microbiol., 10 : 437-441.

Pasculle, A.W., Feeley, J.C., Gibson, R.J., Cordes, L.G., Myerowitz, R.L., Patton, C.M., Gorman, G.W., Carmack, C.L., Ezzell, J.W. and Dowling, J.N. 1980. Pittsburgh pneumonia agent : direct isolation from human lung tissue. J. Infect. Dis., 141 : 727-732.

Edelstein, P.H. 1981. Improved semiselective medium for isolation of *Legionella pneumophila* from contaminated clinical and environmental specimens. J. Clin. Microbio., 14 : 298-303.

Wadowsky, R.W. and Yee, R.B. 1981. Glycine-containing selective medium for isolation of *Legionellaceae* from environmental specimens. Appl. and Environ. Microbiol., 42 : 768-772.

Dennis, P.J.L., Bartlett, C.L.R. and Wright, A.E. 1984. Comparison of isolation methods for *Legionella* spp. In Thornsby, C. et al. (ed.) Legionella : Proceedings of the 2nd International symposium Washington D.C. Am. Soc. Microbiol., p 294-296.

ISO 11731. 1998. Qualité de l'eau – Recherche et dénombrement des *Legionella*.

NF EN ISO 11731-2 (T 90-430). Juillet 2008. Qualité de l'eau. Recherche et dénombrement des *Legionella*. Partie 2 : Méthode par filtration directe sur membrane pour les eaux à faible teneur en bactéries.

11 OTHER INFORMATIONS

The information provided on the package or instructions takes precedence over the formulations or instructions described in this document and are susceptible to modification at any time, without warning.

Code document : BM072/A/2012-05 : 1.
Creation data : 05-2012




CONFIRMATION AGARS FOR *LEGIONELLA* BCYE α AND BCYE α WITHOUT CYSTEINE

Methodology:

- Re-streak 3 characteristic *Legionella* colonies from plate of GVPC agar onto BCYE α (BM072) agar and BCYE α agar without cysteine (BM073).
- Incubate at 36 \pm 1°C for 2 days.

Results:

Growth on BCYE α agar and BCYE α agar without cysteine after incubate at 36°C for 2 days



Legionella pneumophila on BCYE α agar.

Microorganisms	Characteristic colonies on BCYE α agar	Characteristic colonies on BCYE α without cysteine agar
<i>Legionella</i>	White to gray colonies with a smooth surface ; may present a ground glass appearance under a binocular scope	Absence of colony

Product codes:

BM07208 : BCYE α agar - 20 plates

BM07308 : BCYE α agar without cysteine - 20 plates