NOVOBIOCIN SELECTIVE SUPPLEMENT S

ISOLATION OF *SALMONELLA* AND OF *ESCHERICHIA COLI* O157

1 INTENDED USE

Novobiocin is used as a selective supplement in several applications :

- Muller Kauffmann, MKTTn broths, MSRV medium for the selective isolation of motile *Salmonella*.
- Modified casein-soy enrichment broth (mTSB) for the detection of *Escherichia coli* O157.

Novobiocin is proposed in 2 formats, in freeze-dried, sterile form.

2 TYPICAL COMPOSITION

BS033, for one vial :

- Novobiocin 10,0 mg

BS056, for one vial :

- Novobiocin 40,0 mg

3 RECONSTITUTION

BS033 :

- Rehydrate the pellet by aseptically adding 5 mL of sterile distilled water.
- Turn end over tend to dissolve. Avoid frothing the solution.

BS056 :

- Rehydrate the pellet by aseptically adding 20 mL of sterile distilled water.
- Turn end over tend to dissolve. Avoid frothing the solution.

4 QUALITY CONTROL

Aspect : white pellet, giving a colorless and limpid solution after reconstitution.

Typical culture response in MSRV medium, after 24 hours of incubation at 41,5 °C (NF EN ISO 11133) :

Microorganisms		Growth
<i>Salmonella</i> Typhimurium	WDCM 00031	White opaque culture, ≥ 30 mm
<i>Salmonella</i> Enteritidis	WDCM 00030	White opaque culture, ≥ 30 mm
<i>Salmonella</i> Typhimurium	WDCM 00031	White opaque culture, ≥ 30 mm
+ <i>Escherichia coli</i>	WDCM 00013	
+ <i>Pseudomonas aeruginosa</i>	WDCM 00025	
<i>Escherichia coli</i>	WDCM 00013	Inhibited
<i>Enterococcus faecalis</i>	WDCM 00087	Inhibited

5 STORAGE / SHELF LIFE

Freeze-dried Supplements : 2-8 °C.

The expiration date is indicated on the label.

Rehydrated Supplements (*) : 30 days at 2-8 °C.

(*) Benchmark value determined under standard preparation conditions, following manufacturer's instructions.

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8 ADDITIONAL INFORMATION

The information provided on the labels take precedence over the formulations or instructions described in this document and are susceptible to modification at any time, without warning.

Document code : NOVOBIOCINSUPPLEMENTS_ENv7
Creation date : 06-2003
Updated : 05-2016
Origin of revision : General update.