

# RHAPSODY<sup>®</sup>

ENUMERATION OF *PSEUDOMONAS* SPP. IN HUMAN FOOD PRODUCTS  
AND ENVIRONMENTAL SAMPLES OF PRODUCTION AREA

## RELIABLE

Validated method by AFNOR Certification according to NF EN ISO 16140

## PERFORMANCE

Detection of all *Pseudomonas* spp. and total inhibition of secondary flora

## EASY

The blue to blue-green colonies of *Pseudomonas* spp. are easily identifiable

## ECONOMIC

Direct reading without confirmation

## RAPID

Detection and enumeration in only 48 hours

# RHAPSODY

RHAPSODY Agar® allows the detection and the enumeration of *Pseudomonas* spp. in human food products and environmental samples of production area




Validated by AFNOR Certification under the ref. BKR 23/09-05/15 A (meat products) and BKR 23/09-05/15 B (dairy products)



(x) g of sample in 9 (x) mL of diluent <sup>1</sup>

0.1 mL on RHAPSODY Agar®  
(spreading or spiral method) <sup>2, 3</sup>



⊕ Incubation

48 ± 3 h  
30 ± 1 °C


**D +2**

Reading<sup>4</sup> **WITHOUT CONFIRMATION**

Enumeration of **BLUE to BLUE-GREEN**  
colonies of *Pseudomonas* spp.

## To know

With **RHAPSODY Agar®** get result in 2 days and eliminate the confirmation tests specified on the standards.



<sup>1</sup> BPW, Tryptone-salt or any other diluent recommended by the corresponding part of NF EN ISO 6887 standard.

<sup>2</sup> The enumeration limit can be reduced by a factor of 10 by inoculating 1 mL onto the surface of 3 Petri dishes of 90 mm diameter.

<sup>3</sup> The membrane filtration method may be used for environmental samples (out of validation fields).

<sup>4</sup> Reading can be realized after 45 to 72 hours of incubation.

Please refer to the technical data sheet for more information.

## To order

**RHAPSODY Agar® pre-poured**  
BM16708 – 20 Petri dishes (Ø90 mm)

**RHAPSODY Agar® dehydrated**  
BK214HA – 500 g vial  
BS08908 – Suppl. 10 vials qsf 500 mL