




	Dilumat 4	Dilumat S	HiVDS
			
	Great traceability	Great value	High Volume
Dilution Speed	Standard: 1 liter / min With Turbo pump: 1,875 ml / min	750 ml / min	5.8 L / min
Standard pump(s)	1	1	1 double
Additional pump(s)	1	1	0
Dilution Factor	1/2 to 1/50	1/2 to 1/99	Distribution only
Weighing accuracy	< .1g of sample weight	< 1% of sample weight from 5 g < 5% of total weight starting at 10 g	0.50%
Minimum sample weight	1 g	1 g	1 g
Maximum sample weight	1.5 Kg	2.5 Kg	5 Kg
Dispensing mechanism	Automatic arm	Automatic rotating arm	Manual
Bar Code Reader	Optional	No	No
Traceability	Alphanumerical Printer + PC output traceability Diluent	Numerical Printer + PC output	Numerical Printer + PC output
Dimensions L X D X H (mm)	500 x 500 x 500	442 x 450 x 455	442 x 450 x 455
Weight (Kg)	12	11	12
Programming	16	1 per pump	1 per pump
Voltage	115-230 V, 50-60 HZ	115-230 V, 50-60 HZ	115-230 V, 50-60 HZ
Ticket Printer	Optional	Optional	Optional
Bowl stand	Optional	Optional	Standard
Supplier	