






	EASY MIX	SMASHER	PULSIFIER	MINI MIX	JUMBO MIX
					
	Low price	Best value	Specialized use	Smaller samples	Larger samples
Type of blending	Standard	Peristaltic	Ultra sounds	Standard	Standard
Noise level	60 to 70 dBa	48 dBa	60 to 70 dBa	60 to 70 dBa	60 to 70 dBa
Speed	Fixed (240 stroke / min)	3 Speeds (500, 560, 620 stroke / min)	N/A	Variable (360 - 540 strokes / min)	Variable (180 - 420 strokes / min)
Timer	Yes	Yes (10 sec to 3 min)	Yes	Yes (10 to 360 sec.)	Yes (30 to 210 sec.)
Weight	18 Kg	21 Kg	3.5 Kg	9 Kg	45 Kg
Max Blending capacity	80 - 400 ml	80 - 400 ml	80 - 400 ml	100 ml	200 - 2000ml
Waste drawer	Optional	Standard	Standard	Optional	N/A
Heating Capacity	No	No	No	No	Optional
Dimensions (W-H-D)	600 X 400 X 265 mm	400 X 330 X 380 mm	320 X 320 X 460 mm	200 X 200 X 350 mm	450 X 390 X 420 mm