

## BagOpen®

Bag opener/holder

Ref 211 010 (80/100 mL) - Ref 211 040 (400 mL) - Ref 211 350 (3500 mL)



Manufactured under:



BagOpen® allows to easily open the bag and keep it open in vertical position without any contact with the sample.

### TECHNICAL SPECIFICATIONS

- Stable and robust base
- Polycarbonate sides
- Resistant up to 65°C
- Available for 80/100, 400, 3500 mL bags
- Designed and made in France

#### BagOpen® 80/100

Ref: 211 010

- Available for 80- and 100-mL bags
- Dimensions (w x d x h): 9 x 5 x 17 cm
- Net weight : 0.2 kg
- Box (w x d x h): 23 x 12 x 5 cm, gross weight: 0.20 kg

#### BagOpen® 400

Ref: 211 040

- Available for 400- and 2000-mL bags
- Dimensions (w x d x h): 15 x 7.7 x 25.2 cm
- Net weight : 0.4 kg
- Box (w x d x h): 26 x 16 x 6.2 cm, gross weight: 0.54 kg

#### BagOpen® 3500

Ref: 211 350

- Available for 3500 mL bags
- Dimensions (w x d x h): 30 x 18 x 44.1 cm
- Net weight : 2 kg
- Box (w x d x h): 50 x 40 x 40 cm, gross weight: 2.1 kg

### CAUTION

- Do not put in an autoclave