

BagRack®

Storage rack for bags

Ref 221 010 (80/100 mL) - Ref 221 040 (400 mL) - Ref 221 350 (3500 mL - 4 compartments)
Ref 221 352 (3500 mL - 2 compartments)



Manufactured under:



BagRack® enables the clean and methodical storage of bags.

TECHNICAL SPECIFICATIONS

- All stainless steel
- Autoclavable
- Stackable
- Available for 80/100, 400, 2000/3500 mL bags
- Designed and made in France

BagRack® 80/100

Ref: 221 010

- For 11 bags 80/100 mL
- Dimensions (w x d x h): 27.6 x 12.9 x 15.3 cm, net weight: 0.67 kg
- Box (w x d x h): 28.5 x 14.5 x 16 cm, gross weight: 0.7 kg

BagRack® 400

Ref: 221 040

- For 10 bags 400 mL
- Dimensions (w x d x h): 35.5 x 21 x 20 cm, net weight: 1.01 kg
- Box (w x d x h): 37 x 22 x 21 cm, gross weight: 1.37 kg

BagRack® 3500 - 4 compartments

Ref: 221 350

- For 4 bags of 2000/3500 mL
- Dimensions (w x d x h): 43 x 37 x 37 cm, net weight: 2.53 kg
- Box (w x d x h): 45 x 39 x 39 cm, gross weight: 3.44 kg

BagRack® 3500 - 2 compartments

Ref: 221 352

- For 2 bags of 2000/3500 mL
- Dimensions (w x d x h): 35.4 x 21 x 20 cm, net weight: 1.11 kg
- Box (w x d x h): 37 x 23 x 22 cm, gross weight: 1.54 kg

OPTIONS / ACCESSORIES

- Ref. 221 042 DripTray®: Stainless steel drip tray for BagRack® 400 and BagRack® 3500 – 2 compartments
Box (w x d x h): 36 x 22.2 x 3 cm, gross weight: 1 kg