


15mL Polypropylene Conical Centrifuge Tubes with Printed Graduations

Item# 6285/6286/6284



Item# 6285

Description

General purpose centrifuge tube, 500/case. 15mL capacity. Blue screwcap. Black printed graduations. White printed writing space. Available in sterile with attached caps and non-sterile with separate caps.

Features

- Tube molded from chemical resistant, hydrophobic polypropylene (PP)
- Cap molded from polypropylene (PP)
- Conical style bottom
- Printed black graduations (0.5mL graduations)
- White writing area

Products

Item#	Description	Packaging
6285	15mL, PP, sterile, 25 per bag	500 per case
6286	15mL, PP, sterile, 50 per rack	500 per case
6284	15mL, PP, non-sterile, bulk, with separate blue screwcap	500 per case



Item# 6286 in 50 place rack

15mL Polypropylene Conical Centrifuge Tubes with Printed Graduations


Item# 6285/6286/6284

Technical Specifications

Material: Tube: Hydrophobic Polypropylene (PP)
Cap: Polypropylene (PP)
Sterilization: Ethylene Oxide (EO)
Autoclavable: Yes (without cap)
Temperature range: -20°C to +121°C
Spin rate: 8500 RCF

Product Certifications

- ISO 13485
- CE



Customer service, sales and technical support:

Phone 1 (800) 394-4562 • 1 (201) 599-1400

Fax 1 (201) 599-1406

E-mail mail@globescientific.com

Globe Scientific Inc.

610 Winters Avenue • Paramus, NJ 07652

www.globescientific.com

© Copyright 2021 Globe Scientific Inc. All rights reserved.