

## PRODUCT INFORMATION

Your Smarter Culture Media Choice



### Multi Peptone

Cat . No . M13 - 128

Date of issue  
2/22/16

## DESCRIPTION

Multi Peptone is a general purpose growth peptone used for the cultivation of most bacteria, fungi, and some protozoa.

## STORAGE

Store the sealed bottle containing the dehydrated medium at 2 to 30°C. Once opened and recapped, place the container in a low humidity environment at the same storage temperature. Protect it from moisture and light. The dehydrated medium should be discarded if it is no longer free flowing, or if the color has changed.

3651 Clipper Mill Rd. · Baltimore, MD 21211 · Phone (410) 467-9983 · Fax (410) 467-5088  
[www.alphabiosciences.com](http://www.alphabiosciences.com) · [info@alphabiosciences.com](mailto:info@alphabiosciences.com)