


Choose your colony counter

Scan® range



Scan® 300

Best value for money

Minimum colony size: 0.1 mm

Scan® 500

The reference for daily analyses

Minimum colony size: 0.1 mm

Scan® 1200

High performance levels on a wide range of media

Minimum colony size: 0.05 mm

Scan® 4000

Top-of-the-line technology with HD camera and unique lighting features

Minimum colony size: 0.05 mm

	Reference	436 300	436 000	437 000	438 000
Top features	LED lighting	✓	✓	✓	✓
	HD camera (megapixel)	✓	✓	✓	✓
	High sensitivity sensor (CCD)	-	-	✓	✓
	Ultra HD camera (5 megapixels)	-	-	-	✓
	White LED Dome lighting	-	-	-	✓
Colony counting	Pour and surface mode	✓	✓	✓	✓
	Spiral® mode	✓	✓	✓	✓
	Chromogenic agar	-	✓	✓	✓
	Petrifilm™ Compact Dry™ Sanita-kun™	-	-	✓	✓
	Filtration membrane	-	-	✓	✓
	Round Petri dishes up to ø 150 mm	-	-	-	✓
Inhibition zone	90 mm round Petri dishes	-	✓	✓	✓
	120 mm square Petri dishes	-	-	-	✓