

TECHNICAL DATA SHEET

PENICILLIN AND PIMARICIN AGAR (PPA)

ENUMERATION OF *PSEUDOMONAS* spp

1 INTENDED USE

Penicillin and Pimaricin Agar (PPA) is a selective medium for the enumeration of pigmented and non-pigmented psychrotrophic *Pseudomonas* in milk, milk products and in dairy environment samples.

The typical composition responds to that defined in the standard XP ISO/TS 11059.

2 PRINCIPLES

Pancreatic digest of gelatin and Tryptone are the nutrient substrates required for the rapid multiplication of *Pseudomonas*.

Magnesium chloride and potassium sulfate maintain the osmotic balance of the medium and stimulate the production of pyocyanin (a blue, non-fluorescent pigment soluble in water and in chloroform) of *Pseudomonas*.

The penicillin G and pimaricin concentrations inhibit most contaminating bacteria contaminants.

3 TYPICAL COMPOSITION

The composition can be adjusted in order to obtain optimal performance.

For 1 liter of complete media :

- Pancreatic digest of gelatin	16,0 g
- Tryptone	10,0 g
- Potassium sulfate.....	10,0 g
- Magnesium chloride	1,4 g
- Penicillin G	10 ⁵ IU
- Pimaricin.....	10,0 mg
- Bacteriological agar.....	12,0 g

pH of the ready-to-use media at 25 °C : 7,2 ± 0,2.

4 INSTRUCTIONS FOR USE

- Prepare the initial suspension of the product to be analyzed following the XP ISO/TS 11059 standard instructions.
- To the surface of the pre-poured media plates (BM156), transfer 0,1 mL of the stock solution and successive serial dilutions.
- Spread the inoculum over the surface with the help of a sterile triangle or « hockey stick ».
- Incubate at 25 ± 1 °C for 48 ± 2 hours.

✓ **Inoculation :**
0,1 mL on surface

✓ **Incubation :**
48 ± 2 h at 25 ± 1 °C

5 RESULTS

Only count plates containing no more than 150 colonies maximum.

Randomly pick 5 colonies on each plate to perform their purification on 2% Nutrient Agar (BK185HA or BM11808) before confirmation of the genus *Pseudomonas* with the oxidase and glucose fermentation tests (Glucose agar : BM09908).

Pseudomonas species are characterized by the positive reaction to the oxidase test and the absence of glucose fermentation.

6 QUALITY CONTROL

Complete media : whitish agar.

Typical culture response on complete media after 48 hours of incubation at 25 °C (NF EN ISO 11133) :

Microorganisms	Growth (Productivity Ratio : P_R)
<i>Pseudomonas aeruginosa</i>	WDCM 00025
<i>Pseudomonas fluorescens</i>	WDCM 00115
<i>Escherichia coli</i>	WDCM 00013

7 STORAGE / SHELF LIFE

Pre-poured media in plates : 2-8 °C.

The expiration date is indicated on the label.

8 PACKAGING

Pre-poured media in Petri plates (\varnothing 90 mm) :

20 plates BM15608

9 BIBLIOGRAPHY

XP ISO/TS 11059. Octobre 2009. Lait et produits laitiers. Méthode de dénombrement des *Pseudomonas* spp..

10 ADDITIONAL INFORMATION

The information provided on the labels take precedence over the formulations or instructions described in this document and are susceptible to modification at any time, without warning.

Document code : PPA AGAR_ENv2

Creation date : 12-2010

Updated : 11-2015

Origin of revision : General update.