

Your local distributor

## interscience in the world

### PARIS

Phone: +33 (0)1 34 62 62 61 - Email: [info@interscience.com](mailto:info@interscience.com)

### FRANKFURT

Phone: +49 611 7238 7770 - Email: [sales.germany@interscience.com](mailto:sales.germany@interscience.com)

### BOSTON

Phone: +1 781 937 0007 - Email: [sales.usa@intersciencelab.com](mailto:sales.usa@intersciencelab.com)

### SHANGHAI

Phone: +86 (0)21-64739390 - Email: [sales.china@interscience.cn](mailto:sales.china@interscience.cn)

### SINGAPORE

Phone: +65 6977 7232 - Email: [sales.asia@interscience.com](mailto:sales.asia@interscience.com)

### TOKYO

Phone: +81 3 6712 9715 - Email: [sales.japan@interscience.com](mailto:sales.japan@interscience.com)



Anabac.B. 0424. Pictures and information in this brochure are not contractually binding. Cineq<sup>®</sup>, Cineq<sup>®</sup> Anti-algae are trademarks. This document is the translation of the french brochure and is not contractually binding. INTERSCIENCE reserves the right to change or improve the specifications of its products without notice. Please visit [www.interscience.com](http://www.interscience.com) for current updates and additional information. RCS 950 356 220 Versailles. INTERSCIENCE SARL F78860.

# Anabac<sup>®</sup>

Autoclave deodorant



interscience

# Give some air to your lab!

Laboratories are very concerned about hygiene, cleanliness and comfort of work. Try **Anabac**<sup>®</sup>, autoclave deodorant, for a perfumed touch in your lab. **Anabac**<sup>®</sup> autoclave deodorant bring the discrete final touch, the refinement of cleanliness.

## Classic



Eucalyptus fragrance

Ref 320 100



## Poma



Apple fragrance

Ref 320 200



## Citrus



Lemon fragrance

Ref 320 300




## Peach




Peach fragrance


Ref 320 500



1. During sample sterilization, bad smell can emerge from the autoclave.



2. Place one capsule (or 2 if the odor is strong) in the water at the bottom of the autoclave or in a beaker of water.



3. After a sterilization cycle, the capsule melts, allowing the perfume to be released.

## Advantages

- Work comfort
- Compatible with any autoclave
- No effect on the autoclaving
- Pleasant perfume
- Biodegradable

## Applications

- For any kind of autoclave

## Specifications

- Pot of 100 capsules
- Pot weight: 160 g
- 4 fragrances: - Classic (ref: 320 100)  
- Poma (ref: 320 200)  
- Citrus (ref: 320 300)  
- Peach (ref: 320 500)