

	<b>BagMixer® 400 S</b> 400 mL lab blender	Ref. 025 000
	New generation lab blender: new design and silent blending.  BagMixer® 400 S W is the next generation of lab blender: silent, easy-to-use and robust. With great new features, along with the silent blending, it is adapted to all kinds of applications and with a guarantee of optimal bacterial extraction, it is a great tool for lab analyses avoiding all risks of cross-contamination.	

## ADVANTAGES

- ▶ Super silent homogenizing
- ▶ Quick homogenizing
- ▶ Easy cleaning and user-friendly features

## TECHNICAL SPECIFICATIONS

- ▶ 400 ml lab blender
- ▶ Useful volume: 50 - 400 mL
- ▶ Silent blending: less than 48 db
- ▶ Noise cancelling casing
- ▶ Paulstra™ feet: great stability during blending
- ▶ Adjustable blending speed (Slow - 4 strokes/sec; Medium - 6 strokes/sec; Standard - 8 strokes/sec; Quick - 10 strokes/sec)
- ▶ Adjustable blending time (1 second - 1 hour or open-running)
- ▶ Adjustable blending power: soft, standard, strong with 28 kg pressure maximum (Patented)
- ▶ Blending Light Code: green for Ready, orange for Blending, red for Leak detected
- ▶ EnergySaver: 50% less energy consumed with Brushless motor
- ▶ Auto Sleep function: switch to sleeping mode after 30 mn of inactivity
- ▶ Total Access system: total opening of the door for full access of the chamber (Patented)
- ▶ Lighted polycarbonate blending chamber with round corners: easy cleaning
- ▶ Side-by-side paddle stop
- ▶ Easy Cleaning with Click & Clean® removable paddles (Patented)
- ▶ Integrated security drip tray
- ▶ Ambidextrous handle
- ▶ Multi-function digital screen
- ▶ Q-Tight closing: pressure to avoid bag leaks from top
- ▶ Waterproof casing: all stainless steel 304L
- ▶ Motor: heavy-duty brushless motor
- ▶ Power: ~100-240 V / 50-60 Hz
- ▶ Size: 45 x 29 x 37 cm
- ▶ Weight: 26 kg
- ▶ 3 year guarantee (After registration of the Guarantee form)
- ▶ Shock absorbers life-time guarantee
- ▶ Made in France under ISO 9001:2008, CE/RoHS/WEEE
- ▶ In compliance with: ISO 7218
- ▶ Ref. 025 000

---

## DELIVERED WITH

---

- ▶ Set of filter bags
- ▶ Security drip tray
- ▶ Power cord
- ▶ User manual
- ▶ Quick user guide

---

## OPTIONS / ACCESSORIES

---

- ▶ **Ref. 211 040** BagOpen<sup>®</sup> 400, bag opener stand for holding bags
- ▶ **Ref. 221 040** BagRack<sup>®</sup> 400, stainless steel stackable storage racks for 10 bags
- ▶ **Ref. 231 040** BagClip<sup>®</sup> 400, closing clips for lab blender bags 400 mL
- ▶ **Ref. 251 091** BagPipet<sup>®</sup>, automatic pipet, fixed volume 0.1/0.9/1 ml, ejector for BagTips<sup>®</sup>

---

## CONSUMABLES

---

- ▶ **Ref. 111 425** BagFilter<sup>®</sup> 400 P, sterile side filter blender bags
- ▶ **Ref. 122 025** BagPage<sup>®</sup> + 400, sterile full-page filter blender bags
- ▶ **Ref. 132 025** BagLight<sup>®</sup> 400 PolySilk<sup>®</sup>, plain blender bags
- ▶ **Ref. 145 040** RollBag<sup>®</sup> 1300, sampling blender bags
- ▶ **Ref. 252 019** BagTips<sup>®</sup> REGULAR, 19 cm sterile tips for BagPipet<sup>®</sup> or other pipet
- ▶ **Ref. 252 024** BagTips<sup>®</sup> JUMBO, 24 cm sterile tips for BagPipet<sup>®</sup> or other pipet