


|   |   |              |
|---|---|--------------|
|  | <b>BagMixer® 400 SW</b><br>400 mL lab blender   | Ref. 025 100 |
|   | New generation lab blender: silent blending and new design<br>BagMixer® 400 SW is the next generation of lab blenders with a new design: silent, easy-to-use and robust. With great new features, along with the silent blending, it is adapted to all kinds of applications and with a guarantee of optimal bacterial extraction, it is a great tool for lab analyses avoiding all risks of cross-contamination. |              |

## ADVANTAGES

- ▶ Super silent homogenizing
- ▶ Quick homogenizing with window door
- ▶ LiquidSensor (Patented)
- ▶ Easy cleaning and user-friendly features

## TECHNICAL SPECIFICATIONS

- ▶ 400 ml lab blender
- ▶ Useful volume: 50 - 400 mL
- ▶ Silent blending: less than 48 db
- ▶ Noise cancelling casing
- ▶ Double layer polycarbonate window
- ▶ Paulstra™ ventouses feet: great stability during blending
- ▶ Adjustable blending speed (Slow - 4 strokes/sec; Medium - 6 strokes/sec; Standard - 8 strokes/sec; Quick - 10 strokes/sec)
- ▶ Adjustable blending time (1 second - 1 hour or open-running)
- ▶ Adjustable blending power: choose your blending power: Soft/Standard/Strong with up to 28 kg pressure maximum per paddle)
- ▶ Total Access system: total opening of the door for full access of the chamber (Patented)
- ▶ Side-by-side paddle stop
- ▶ LiquidSensor: immediate stop in case of liquid leak (Patented)
- ▶ Easy Cleaning with Click & Clean removable paddles (Patented)
- ▶ Integrated security drip tray
- ▶ Non-marking Paulstra™ feet
- ▶ Ambidextrous handle
- ▶ Multi-function digital screen
- ▶ Blending chamber: water-proof, smooth, LED Lighting, White gloss polycarbonate
- ▶ Blending LightCode: green for Ready, orange for Blending, red for Leak detected
- ▶ Energy saver after 30 minutes: auto-sleep after 30 minutes of inactivity
- ▶ Q-Tight closing: 250 kg pressure to avoid bag leaks
- ▶ Waterproof casing: all stainless steel 304L
- ▶ Motor: heavy-duty 180 Watts brushless motor
- ▶ Power: ~100-240 V / 50-60 Hz
- ▶ Size: 45 x 29 x 37 cm
- ▶ Weight: 26 kg
- ▶ Box: size (W x D x H): 57 x 42 x 42 cm, weight: 29.5 kg
- ▶ 3 year guarantee (After registration of the Guarantee form)
- ▶ Window door life-time guarantee
- ▶ Shock absorbers life-time guarantee
- ▶ Made in France under ISO 9001:2008, CE/RoHS/WEEE
- ▶ In compliance with: ISO 7218

- ▶ Delivered with: power cord, user's manual, security drip tray
  - ▶ Ref. 025 100
- 

---

## DELIVERED WITH

---

- ▶ Security drip tray
  - ▶ Power cord
  - ▶ User manual, IQ/OQ/PQ guidelines
  - ▶ Quick user guide
- 

## OPTIONS / ACCESSORIES

---

- ▶ **Ref. 211 040** BagOpen<sup>®</sup> 400, bag opener stand for holding bags
  - ▶ **Ref. 221 040** BagRack<sup>®</sup> 400, stainless steel stackable storage racks for 10 bags
  - ▶ **Ref. 231 040** BagClip<sup>®</sup> 400, closing clips for lab blender bags 400 mL
  - ▶ **Ref. 251 091** BagPipet<sup>®</sup>, automatic pipet, fixed volume 0.1/0.9/1 ml, ejector for BagTips<sup>®</sup>
- 

## CONSUMABLES

---

- ▶ **Ref. 111 425** BagFilter<sup>®</sup> 400 P, sterile side filter blender bags
- ▶ **Ref. 122 025** BagPage<sup>®</sup> + 400, sterile full-page filter blender bags
- ▶ **Ref. 132 025** BagLight<sup>®</sup> 400 PolySilk<sup>®</sup>, plain blender bags
- ▶ **Ref. 145 040** RollBag<sup>®</sup> 1300, sampling blender bags
- ▶ **Ref. 252 019** BagTips<sup>®</sup> REGULAR, 19 cm sterile tips for BagPipet<sup>®</sup> or other pipet
- ▶ **Ref. 252 024** BagTips<sup>®</sup> JUMBO, 24 cm sterile tips for BagPipet<sup>®</sup> or other pipet