


interscience

30, CH. BOIS ARPENTS – 78860 ST-NOM-LA-BRETECHE – FRANCE – T. +33 (0)1 34 62 62 61 – F. +33 (0)1 34 62 43 03
N° Identification TVA: FR52 950 356 220 – E-mail: info@interscience.fr – Website: www.interscience.fr

	MiniMix® P CC® 100 ml lab blender	Ref. 011 230
	Quick and accurate MiniMix® 100 is designed for the blending of small samples Thanks to the adjustable paddles, blending is more accurate and efficient. Cleaning is easy with Click & Clean® removable paddles system and 270° fully opening door.	

ADVANTAGES

- ▶ Silent and resistant
- ▶ Variable blending power
- ▶ Easy cleaning: door opens at 270°

TECHNICAL SPECIFICATIONS

- ▶ All stainless steel
 - ▶ Plain door
 - ▶ Multifunction digital screen
 - ▶ Variable speed (3, 6, 9, 12 strokes/sec.)
 - ▶ Variable time (1 sec. - 1h or ∞)
 - ▶ Variable power, adjustable paddles (patented)
 - ▶ Q-Tight® closing
 - ▶ Click & Clean® removable paddles (patented)
 - ▶ Fully opening door (270°)
 - ▶ Side-by-side paddle stop
 - ▶ Useful volume: 5 – 80 mL
 - ▶ Also available in CoolMix® version
- ▶ Power: 110-220V~ 50-60Hz
Machine: size (W x D x H): 21 x 37 x 20 cm, weight: 13 kg
Box: size (W x D x H): 48.5 x 28.5 x 30 cm, weight: 15.5 kg
- ▶ 3 year guarantee (after registration)
 - ▶ Shock absorbers life-time guarantee
 - ▶ Life-time service
 - ▶ In compliance with: ISO 7218
 - ▶ Designed and made in France
Manufactured under: ISO 9001 V2008, CE, WEEE 2002/96/EC, RoHS 2002/95/EC standards

DELIVERED WITH

- ▶ Power cord
- ▶ User manual, IQ/OQ/PQ guidelines
- ▶ Quick user guide
- ▶ Free sample pack of bags

OPTIONS / ACCESSORIES

- ▶ **Ref. 016 000** Drip tray MiniMix®
- ▶ **Ref. 211 010** BagOpen® 80/100, bag opener stand for holding bags
- ▶ **Ref. 221 010** BagRack® 80/100, stainless steel stackable storage rack for 11 bags of 80/100 mL
- ▶ **Ref. 231 010** BagClip® 80/100, closing clips for lab blender bags 80/100 mL
- ▶ **Ref. 251 091** BagPipet®, autoclavable pipet with fixed volume: 0.1 mL/0.9 mL/1 mL with BagTips®

CONSUMABLES

- ▶ **Ref. 120 825** BagPage® 80, sterile full-surface filter bag
- ▶ **Ref. 121 025** BagPage® 100, sterile full-surface filter bag
- ▶ **Ref. 131 025** BagLight® 100, sterile plain bag
- ▶ **Ref. 145 010** RollBag® 450, sterile sampling bag
- ▶ **Ref. 252 019** BagTips® REGULAR, 19 cm sterile tips for BagPipet® or other pipet
- ▶ **Ref. 252 024** BagTips® JUMBO, 24 cm sterile tips for BagPipet® or other pipet