interscience

Scan® 50 pro

Manual colony counter

Ref. 435 055







FDA BAM ISO

Traçabilité et Haute technologie :







Fabrication:









Scan® 50 pro is a manual colony counter with multiple ergonomic supports (posture, vision, hearing). The Dark Field 2 technology offers an indirect, precise, and contrasted LED lighting of the colonies. The export of the counting result, by USB link, to a computer, ensures the traceability of the data, without preliminary processing by software.

TECHNICAL SPECIFICATIONS

- · Stainless steel frame, very robust
- Multifunction LCD screen
- · Storage hatch for counting grids and white backgrounds
- · Inclined and padded counting surface
- Dark Field 2 technology for high contrast colony lighting
- Adjustable LED intensity: 5 lighting levels and can be turned off Complies with IEC 62471:2006 standard (group 2 risk)
- Counting: 0 to 1 999 UFC
- For all types of support: Petri dishes: Ø 55 90 mm / PetriFilm™, Compact Dry™, MC-Media Pads™, easyPlate™ / filtration membranes
- Touch surface with adjustable sensitivity: 5 levels of detection of a pressure onto support
- Suitable for all types of markers
- 180°-HandPad for left or right forearm
- · Counting tone: choice of 9 different sounds / 5 volume levels
- · Headset for sound reception and isolation of the work environment
- · Counting light signal: green flash emitted by the set of 3 internal LEDs
- · Magnifying glass for colony size (optional)
- Counting keys duplicated on each side of the keyboard, for left or right-handed operators
- USB port: transfer of counting data to computer
- Data export to LIMS, under CSV format (Excel[™]) or text format (Word[™])
- Supply voltage and frequency: 100-240V ~ 50-60Hz
- Maximum power: 20 W
- Fuse: T3,15A L 250V 5x20 mm
- Machine protection index: IP20
- Size (I x p x h): 27.3 x 35.2 x 15.1 cm, weight: 4.45 kg
- Packing (I x p x h): 35 x 45 x 40 cm, weight : 7.05 kg
- · Warranty: 3 years, after registration of the warranty card
- Manufacturing CE, RoHS 2002/95/EC, WEEE 2002/95/EC
- Designed and manufactured in France
- · Availability of spare parts: 10 years
- In compliance with ISO 7218 and FDA BAM standards (Bacteriological Analytical Manual)

1/2

Scan® 50 pro_FT 11/22 Pictures and information are not contractually binding. Scan® is trademark. INTERSCIENCE reserves the right to change or improve the specifications of its products without notice. Please visit www.interscience.com for current updates and additional information. RCS 950 356 220 Versailles. INTERSCIENCE SARL F78860

interscience

SUPPLIED WITH

- 1 power cord
- 1 USB cable
- 1 black marker pen
- 1 white diffusing background
- 1 user manual
- 1 headset

OPTIONS / ACCESSORIES

0. 11011077100200011120	
 Ref. 435 001 	Magnifier: Magnifier + 1 flexible arm for Scan® 50 and Scan® 50 pro
• Ref. 435 025	Wolffhuegel counting grid: 1 transparent and 1 white diffusing backgrounds for Scan® 50 and Scan® 50 pro Ø 90 mm
• Ref. 435 036	Spiral® counting grid (90 mm): easySpiral® counting grid for Scan® 50 and Scan® 50 pro - 1 transparent and 1 white diffusing backgrounds Ø 90 mm
• Ref. 435 038	Circle counting grid: easySpiral® Dilute circle counting grid for Scan® 50 and Scan® 50 pro - 1 transparent and 1 white diffusing backgrounds Ø 90 mm
• Ref. 435 051	Adaptor for Petrifilm™, MC-Media Pads™, Compact Dry™: Adaptor for Petrifilm™, MC-Media Pads™, Compact Dry™ for Scan® 50 and Scan® 50 pro
 Ref. 435 052 	Headset
 Ref. 435 702 	Optical cleaner + wipes
 Ref. 522 000 	Barcode reader

info@interscience.com | sales.germany@interscience.com | sales.usa@interscience.com | sales.usa@interscience.com | sales.asia@interscience.com | sales.asia@